



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name: Kurklatex
Item code: KU
Barcode: N/A
CE Marking: N/A
MDR Risk Classification: Riskclass I

Product summary: Kurklatex
Size / Contents product: 92 x 46 cm
Packaging size: plate, cutted shapes
Packaging unit: per plate
Color: beige/brown
Thickness: 1, 1.5, 2, 3, 4, 5, 6, 8 en 10mm
Shorevalue: 60-65 ShoreA
Implementation: N/A
Storage advise: not in direct sunlight
Maintenance advise: N/A
Sterilisation advise: N/A
Manufacturing date: N/A
Shelf life: N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands
Telefoon +31 (0)416 37 69 87 • **Fax** +31(0)416 37 56 02 • **E-mail** info@medical-leather.nl •
Website www.medical-leather.nl **Rabobank** 1308.95.377 • **BIC Code** RABONL2U • **IBAN nr.** NL63 RABO 0130 8953 77 •
KVK 18124129 • **BTW nr.** NL8145 15 745 B01



medical
leather®

a healthy focus on feet

DECLARATION OF CONFORMITY

The undersigned
C.J. Maas – Leder LEFA BV. / Medical Leather
In reference
To the REACH Regulation 1907/2006/CE
Concerning the Registration, Evaluation, Authorization
And the restriction of chemical substances

DECLARES THAT

- Our company purchases and sells foam sheets through the transformation of polymers, mineral fillers, pigments and chemicals. Therefore, as a user, our company is not bound to any registration
- According to the REACH Regulation polymers are exempted from registration (article 2) as well as most of the raw materials used.
- Through the raw material suppliers and sheet suppliers, we were able to ascertain compliance with the REACH legislation. In fact, the chemical manufacturers, where provided, have taken steps as for regular registration.
 - With regard to the presence of hazardous substances (SVHC) mentioned in **the last list published by ECHA on 25/06/2020** we declare that our products, after curing, as supplied by us, has no hint of SVHC in a concentration above the threshold limit of 0.1%..

Pieve del Cairo, 22/07/2020

REACH Responsible
Pieter Maas

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 • KVK 18124129 • BTW nr. NL8145 15 745 B01

Colour	Natural
Cork granule size	1-3 mm
Binder	Synthetic Rubber

TEST METHOD	PROPERTY	UNIT	VALUE
ASTM F1315	Density	kg/m ³	480
ASTM D 2240	Hardness	shore A	45-65
ASTM F36	Compressibility 400 psi	%	40-60
ASTM F36	Recovery 400 psi (min)	%	75
ASTM F152	Tensile strength (min)	MPa	0.69
ASTM F147	Flexibility		
	Original (F=5)	-	No cracks
	Oven aged, 70h at 100°C (F=16)	-	No cracks
	ASTM n°1 Oil, 70h at 100°C (F=16)	-	No cracks
ASTM D471	Volume change after immersion		
	ASTM n°1 Oil, 70h at 100°C	%	-10 to+10
	ASTM n°3 Oil, 70h at 100°C	%	0 to+40
	ASTM fuel A, 22h at RT	%	0 to+15

Dimensional tolerance:

Length	Width
±2%	±2%

Thickness (mm)	Tolerance (mm)
≤ 1,6	±0,25

General Information:

We recommend storing the product at these conditions: 20±3°C / 60±5%HR