



medical<sup>®</sup>  
leather

a healthy focus on feet

## Product Specific Information Document

<b>Product name:</b>	<b>Lena</b>
Item code:	LE 09
Barcode:	N/A
CE Marking:	N/A
MDR Risk Classification:	Riskclass I
Product summary:	Lena 1,0 – 1,2mm
Size / Contents product:	45-50 sqft
Packaging size:	Hide
Packaging unit:	per hide
Color:	oker, white, grey, red, sky, creme, lime, beige, pink, maron, marine, brown, testa, taupe, blue, yellow, royal, purple, black
Thickness:	1,0 – 1,2mm
Shorevalue:	N/A
Implementation:	non-perforated
Storage advise:	not in direct sunlight
Maintenance advise:	This leather can be cleaned with a mild, non alcoholic, cleaning substance
Sterilisation advise:	N/A
Manufacturing date:	N/A
Shelf life:	N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands  
Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •  
Website [www.medical-leather.nl](http://www.medical-leather.nl) Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 •  
KVK 18124129 • BTW nr. NL8145 15 745 B01



medical  
leather®

a healthy focus on feet

## DECLARATION OF CONFORMITY

The undersigned  
C.J. Maas – Leder LEFA BV. / Medical Leather  
In reference  
To the REACH Regulation 1907/2006/CE  
Concerning the Registration, Evaluation, Authorization  
And the restriction of chemical substances

### DECLARES THAT

- Our company purchases and sells foam sheets through the transformation of polymers, mineral fillers, pigments and chemicals. Therefore, as a user, our company is not bound to any registration
- According to the REACH Regulation polymers are exempted from registration (article 2) as well as most of the raw materials used.
- Through the raw material suppliers and sheet suppliers, we were able to ascertain compliance with the REACH legislation. In fact, the chemical manufacturers, where provided, have taken steps as for regular registration.
  - With regard to the presence of hazardous substances (SVHC) mentioned in **the last list published by ECHA on 25/06/2020** we declare that our products, after curing, as supplied by us, has no hint of SVHC in a concentration above the threshold limit of 0.1%..

Pieve del Cairo, 22/07/2020

REACH Responsible  
Pieter Maas

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail [info@medical-leather.nl](mailto:info@medical-leather.nl) •

Website [www.medical-leather.nl](http://www.medical-leather.nl) Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 • KVK 18124129 • BTW nr. NL8145 15 745 B01

# TECHNICAL SPECIFICATION SHEET

ARTICLE NAME:

**LENA**

RAW MATERIAL: BOVINE LEATHER AUSTRALIA/NEW ZEALAND/EUROPE

AVERAGE SIZE: 4.20 m<sup>2</sup>

THICKNESS: ca. 0.9 mm

Test	Test method	Result
1. Colour fastness to rubbing: Nr. 500 dry rubs Nr. 150 wet rubs Nr. 80 alkaline ph 8	IUF 450 / UNI EN ISO 11640 IUF 450 / UNI EN ISO 11640 IUF 450 / UNI EN ISO 11640	≥ 4 grey scale ≥ 4 grey scale ≥ 4 grey scale
2. Flexing endurance	IUP 20 UNI EN ISO 5402	≥ 30.000
3. Tear strength	IUP 40 UNI EN ISO 3377-1	≥ 20 N
4. Adhesion	IUF 470 UNI EN ISO 11644	≥ 3.50 N/cm
5. Light resistance (xenon test)	IUF 402 UNI EN ISO 105-B02	5 blue scale
6. Fire resistance	UNI EN ISO 1021-1/2 BS 5852 1990 Sec. 4 NI/0-1 California T.B. 117	pass pass pass

All our leathers comply with the strict regulations and EU Directives covering the use and presence of PCP, CFC, Chromo VI, Dimethylfumarate (DMF) and Azo dyestuffs during the tanning process. Due to the natural properties of leather colour variations between one dye lot and another may occur. This technical specification sheet replaces any copies in your possession. The values indicated in this sheet refer to the leather at the time of supply. Once the leather has been processed, cut, marked we can no longer guarantee the same values.

Medical Leather  
Zanddonkweg 6  
5144 NX Waalwijk  
The Netherlands