



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name: Podolux
Item code: EV
Barcode: N/A
CE Marking: N/A
MDR Risk Classification: Riskclass I

Product summary: Podolux
Size / Contents product: 95 x 85 cm, 11,5 cm x 25 m
Packaging size: sheet, cutted shapes, wheel
Packaging unit: per sheet, cutted shapes, wheel
Color: swirl green/purple, swirl pink/purple,
Thickness: 1,2 mm, 1,6 mm
Shorevalue: 30 ShoreA
Implementation: non-perforated
Storage advise: not in direct sunlight
Maintenance advise: N/A
Sterilisation advise: N/A
Manufacturing date: N/A
Shelf life: N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 •

KVK 18124129 • BTW nr. NL8145 15 745 B01



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name: Podolux
Item code: EV
Barcode: N/A
CE Marking: N/A
MDR Risk Classification: Riskclass I

Product summary: Podolux
Size / Contents product: 95 x 85 cm, 11,5 cm x 25 m
Packaging size: sheet, cutted shapes, wheel
Packaging unit: per sheet, cutted shapes, wheel
Color: army, zebra black/wit, zebra red/white, zebra lime/blue, zebra pink/pink
Thickness: 1,2 mm
Shorevalue: 30 ShoreA
Implementation: non-perforated
Storage advise: not in direct sunlight
Maintenance advise: N/A
Sterilisation advise: N/A
Manufacturing date: N/A
Shelf life: N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 •

KVK 18124129 • BTW nr. NL8145 15 745 B01



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name:	Podolux
Item code:	EV
Barcode:	N/A
CE Marking:	N/A
MDR Risk Classification:	Riskclass I
Product summary:	Podolux
Size / Contents product:	95 x 85 cm
Packaging size:	sheet, cutted shapes
Packaging unit:	per sheet, cutted shapes
Color:	water
Thickness:	2 mm
Shorevalue:	30 ShoreA
Implementation:	non-perforated
Storage advise:	not in direct sunlight
Maintenance advise:	N/A
Sterilisation advise:	N/A
Manufacturing date:	N/A
Shelf life:	N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands
Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •
Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 •
KVK 18124129 • BTW nr. NL8145 15 745 B01



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name:	Podolux
Item code:	EV
Barcode:	N/A
CE Marking:	N/A
MDR Risk Classification:	Riskclass I
Product summary:	Podolux
Size / Contents product:	95 x 85 cm
Packaging size:	sheet, cutted shapes
Packaging unit:	per sheet, cutted shapes
Color:	aqua blue and aqua black
Thickness:	1,2 mm
Shorevalue:	30 ShoreA
Implementation:	non-perforated
Storage advise:	not in direct sunlight
Maintenance advise:	N/A
Sterilisation advise:	N/A
Manufacturing date:	N/A
Shelf life:	N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands
Telefoon +31 (0)416 37 69 87 • **Fax** +31(0)416 37 56 02 • **E-mail** info@medical-leather.nl •
Website www.medical-leather.nl **Rabobank** 1308.95.377 • **BIC Code** RABONL2U • **IBAN nr.** NL63 RABO 0130 8953 77 •
KVK 18124129 • **BTW nr.** NL8145 15 745 B01



medical[®]
leather

a healthy focus on feet

DECLARATION OF CONFORMITY

The undersigned

C.J. Maas – Leder LEFA BV. / Medical Leather

In reference

To the REACH Regulation 1907/2006/CE

Concerning the Registration, Evaluation, Authorization

And the restriction of chemical substances

DECLARES THAT

- Our company purchases and sells foam sheets through the transformation of polymers, mineral fillers, pigments and chemicals. Therefore, as a user, our company is not bound to any registration
- According to the REACH Regulation polymers are exempted from registration (article 2) as well as most of the raw materials used.
- Through the raw material suppliers and sheet suppliers, we were able to ascertain compliance with the REACH legislation. In fact, the chemical manufacturers, where provided, have taken steps as for regular registration.
- With regard to the presence of hazardous substances (SVHC) mentioned in **the last list published by ECHA on 25/06/2020** we declare that our products, after curing, as supplied by us, has no hint of SVHC in a concentration above the threshold limit of 0.1%..

Pieve del Cairo, 22/07/2020

REACH Responsible

Pieter Maas

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 •

KVK 18124129 • BTW nr. NL8145 15 745 B01



date last revision 28/02/2020

SAFETY DATA SHEET

1- PRODUCT AND COMPANY IDENTIFICATION DETAILS

1.1 Product name: **ML-CAD 20-30 shore**

1.2 PRODUCER: CMP SRL- Via Matteotti 52-27037 Pieve Del Cairo (PV)

1.1 Company phone number for urgent calls: 0384-831772

2- INGREDIENT/COMPOSITION INFORMATION

SUBSTANCE	N.CAS	
E.V.A. (ethylene/vinyl acetate)	24937-78-8	
Calcium carbonate	471-34-1	
Silicon dioxide	7631-86-9	
Stearic acid	67701-03-5	
Zinc oxide	1314-13-2	
Peroxide	25155-25-3	
Azodicarbonamide	123-77-3	
Zinc stearate	91051-01-3	
For the colors, they are used various combinations of the following materials:		
Organic pigments	//	
Yellow iron oxide	51274-00-1	
Red iron oxide	1309-37-1	
Brown iron oxide	//	
Black iron oxide	1317-61-9	
Bitumen	//	
Titanium dioxide	13463-67-7	
Carbon black	1333-86-4	

3. DANGER IDENTIFICATION

3.1 Main health risks:

No specific danger has been found in normal use

3.2 Main risks for the environment

Not to be disposed in the environment



medical
leather®

a healthy focus on feet

4. FIRST- AIO MEASURES

- 4.1 Skin contact: no particular procedure

5. FIRE PREVENTION MEASURES

- 5.1 Recommended extinguishers : Atom ised water , C02, Plastic foam, Chemical powder, Halons.
- 5.2 Forbidden extinguishers : None
- 5.3 Combustion risks: The vapours produced by thermal decomposition may irritate the skin and mucal areas. Avoid breathing the fumes .
- 5.4 Protection means: Protective clothing complete with autorespiration system

6. MEASURES INTHE EVEN OFACCIDENTAL LEAKAG

Not verifiable

7. HANDLING ANO STORAG

- 7.1 Handling: Handling of the material at room temperature is risk-free:
- 7.2 Storage: away from sunlight , in well-aired places that are not excessively hot
- 7.3 Precautions : avoid open flames, place the equipment on the floor to avoid accumulating static electricity .

8. EXPOSURE /INDIVIDUAL PROTECTION CHECKS

- 8.1 Precautions: Adequatel ventilate the premises in which the product is stored.
- 8.2 Respiratory protection: Not necessary
- 8.3 Eye protection: Not necessary
- 8.4 Exposure limits of the substances making up the sam : none



9. CHEMICAL AND PHYSICAL PROPERTIES

- 9.1 Physical state: solid- Available in sheets from 60x120 cm to 120x190 cm, thickness varies from 2 to 40 mm
- 9.2 Colour: available in various colours
- 9.3 Odour: n. a.
- 9.4 pH: n. a.
- 9.5 Density: from 0.15 g/cm³ to 1.20 g/cm³ based on the articles
- 9.6 Inflammability : n. a.

10. STABILITY AND REACTIVITY

- 10.1 The product is stable and chemically inert at room temperature.

11. TOXICOLOGICAL INFORMATION

- 11.1 This product is not toxic. However, in the event of fire, the fumes that may develop may irritate the respiratory system .

12. ECOLOGICAL INFORMATION

- 12.1 The product is not biodegradable . It does not contain additives based on lead, mercury, cadmium and hexavalent chrome. It does not present risks for the water table.
USE THE PRODUCT ACCORDING TO GOOD WORKING PRACTICE AND DISPOSE OF THE PRODUCT IN AN ENVIRONMENTALLY FRIENDLY WAY AFTER USE.

13. DISPOSAL CONSIDERATIONS

- 13.1 Proceed in accordance with current local and national regulations.



14. TRANSPORT INFORMATION

14.1 This product is not considered dangerous goods.

15. REGULATORY INFORMATION

- 15.1 In accordance with D.M. 28/1/92 L-D, this product is not considered dangerous goods.
- 15.2 Where applicable, reference should be made to the following regulations:
D.P.R. 303/56 (sanitary checks) - circulars no. 46/79 and 61/81 - Ministry of Labour.
D.P.R. 175/88 (Seveso Directive), Enclosure 2°, 3° and 4°.

16. FURTHER INFORMATION

- 16.1 Main bibliographic sources: safety sheet for the substances used in the manufacture of this product
The above information is based on our best information to our knowledge on the date indicated above. It refers solely to the product indicated and does not constitute a guarantee of any special qualities. The user must make sure of that the information is suitable and complete in relation to the intended use.
This fact sheet annuls and replaces all previous editions.