



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name:	Sanino
Item code:	SA N
Barcode:	N/A
CE Marking:	N/A
MDR Risk Classification:	Riskclass I
Product summary:	Sanino
Size / Contents product:	150 x 100 cm
Packaging size:	sheet
Packaging unit:	per sheet
Color:	anthracite, beige, brown
Thickness:	approx. 1,5mm
Shorevalue:	N/A
Implementation:	N/A
Storage advise:	not in direct sunlight
Maintenance advise:	N/A
Sterilisation advise:	N/A
Manufacturing date:	N/A
Shelf life:	N/A

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands
Telefoon +31 (0)416 37 69 87 • **Fax** +31(0)416 37 56 02 • **E-mail** info@medical-leather.nl •
Website www.medical-leather.nl **Rabobank** 1308.95.377 • **BIC Code** RABONL2U • **IBAN nr.** NL63 RABO 0130 8953 77 •
KVK 18124129 • **BTW nr.** NL8145 15 745 B01



medical
leather®

a healthy focus on feet

DECLARATION OF CONFORMITY

The undersigned
C.J. Maas – Leder LEFA BV. / Medical Leather
In reference
To the REACH Regulation 1907/2006/CE
Concerning the Registration, Evaluation, Authorization
And the restriction of chemical substances

DECLARES THAT

- Our company purchases and sells foam sheets through the transformation of polymers, mineral fillers, pigments and chemicals. Therefore, as a user, our company is not bound to any registration
- According to the REACH Regulation polymers are exempted from registration (article 2) as well as most of the raw materials used.
- Through the raw material suppliers and sheet suppliers, we were able to ascertain compliance with the REACH legislation. In fact, the chemical manufacturers, where provided, have taken steps as for regular registration.
 - With regard to the presence of hazardous substances (SVHC) mentioned in **the last list published by ECHA on 25/06/2020** we declare that our products, after curing, as supplied by us, has no hint of SVHC in a concentration above the threshold limit of 0.1%..

Pieve del Cairo, 22/07/2020

REACH Responsible
Pieter Maas

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 • KVK 18124129 • BTW nr. NL8145 15 745 B01

SANINO

Composition		100% PU 70% PES 30% CO
Width of roll		140 cm +/- 2 cm
Weight		650 g/lm +/- 10%
Abrasion resistance	EN ISO 5470-2	> 100.000 cycles
Thickness		1,00 mm +/- 0,03 mm
Rubbing solvents	ISO 105-X12:2001	5
Rubbing fastness in humid state	ISO 105-X12:2001	5
Color fastness to light	ISO 105 B02	4-5
Resistance to seam slippage	EN ISO 13936-2	2,5 mm warp 3 mm weft
Tear growth strenght		
Warp	EN ISO 4674-1	105,8 N
Weft		73,2 N
Tensile strenght		
Warp	EN ISO 1421-1	908 N
Weft		593 N
Fire resistance		on request

CLEANING



1. The results of a test were collected by various samples and therefore presenting average indicative values.
2. Colours of fabrics may very slightly away from samples placed in printed and digital catalogs or placed in other media.
3. Color shade may slightly vary from dye-lot to dye-lot.
4. Penetration of dyes fom fabrics with dark and saturated colors, to various bright surfaces is inevitable.
5. Pilling effect is a natural feature of fabrics, completely avoided is not possible.
6. OCTO company don't take any responsibility of any damages resulting fom improper use and cleaning.
7. Mechanical damage is not subject to complaint.