



medical[®]
leather

a healthy focus on feet

Product Specific Information Document

Product name:	SMF
Item code:	SM
Barcode:	N.A.
CE Marking:	N.A.
MDR Risk Classification:	Riskclass I
Product summary:	SMF
Size / Contents product:	137 x 72 cm
Packaging size:	per sheet
Packaging unit:	Per sheet
Color:	taupe / black
Thickness:	0,8 mm
Shorevalue:	N/A
Implementation:	perforated
Storage advise:	not in direct sunlight
Maintenance advise:	N.A.
Sterilisation advise:	N.A.
Manufacturing date:	N/A
Shelf life:	N.A.

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 •

KVK 18124129 • BTW nr. NL8145 15 745 B01



DECLARATION OF CONFORMITY

The undersigned

C.J. Maas – Leder LEFA BV. / Medical Leather

In reference

To the REACH Regulation 1907/2006/CE

Concerning the Registration, Evaluation, Authorization

And the restriction of chemical substances

DECLARES THAT

- Our company purchases and sells foam sheets through the transformation of polymers, mineral fillers, pigments and chemicals. Therefore, as a user, our company is not bound to any registration
- According to the REACH Regulation polymers are exempted from registration (article 2) as well as most of the raw materials used.
- Through the raw material suppliers and sheet suppliers, we were able to ascertain compliance with the REACH legislation. In fact, the chemical manufacturers, where provided, have taken steps as for regular registration.
- With regard to the presence of hazardous substances (SVHC) mentioned in **the last list published by ECHA on 25/06/2020** we declare that our products, after curing, as supplied by us, has no hint of SVHC in a concentration above the threshold limit of 0.1%..

Pieve del Cairo, 22/07/2020

REACH Responsible

Pieter Maas

Postbus 95, 5140 AB Waalwijk • Zanddonkweg 6, 5144 NX Waalwijk • The Netherlands

Telefoon +31 (0)416 37 69 87 • Fax +31(0)416 37 56 02 • E-mail info@medical-leather.nl •

Website www.medical-leather.nl Rabobank 1308.95.377 • BIC Code RABONL2U • IBAN nr. NL63 RABO 0130 8953 77 • KVK 18124129 • BTW nr. NL8145 15 745 B01

Op al onze overeenkomsten zijn van toepassing de leverings- en betalingsvoorwaarden voor leder en andere artikelen, gedeponeerd bij de Kamer van Koophandel te Tilburg.

SMF

Weight (gr/square mt) :	350±25	Thickness (mm) :	0,70÷0,80	Width (cm) :	
Abrasion resistance : 102.400 cycles (Martindale BS 5690)					
Hot test: 160 °C 30 seconds					
Flexion resistance : 20 °C >200000 cycles (Uni 4818/13) -10 °C >100000 cycles -20 °C 50000 cycles					

		Warp	Weft
Tensile strenght : (UNI 4818/6)	N/5cms	500±50	400±50
Extention at break (UNI 4818/6)	%	60±6	90±9
Tear resistance: (UNI 4818/9)	N	30±3	35±3
Stitching : (UNI 4818/11)	N/5cms	350±35	300±30
Skin adhesion : (UNI 4818/10)	N/cm	12±2	12±2

Injectable with :			
Leakproofness test	(UNI 5122)	hydrostatic pressure limit h.=mt.2 (evaluation) :	Sufficiente
Surface resistance to cold liquids UNI EN 12720-UNI 10944			

Composition % :	support	46 %	100%PA
	Resin	54 %	PU

Average datas with usual tollerances. They could be modified without warning.